

6-Axis-Parallel Kinematic Positioning Systems

HXP 100 Hexapods



The HXP100-MECA 6-Axis Hexapod is a parallel kinematic motion device that provides six degrees of freedom: X, Y, Z, pitch, roll, and yaw. The HXP100 has long travel capability and is an effective solution to complex motion applications that demand high load capacity and accuracy in up to six independent axes. Contrary to the image of being complex and highly priced, the HXP100 hexapod is not only affordable but extremely easy to use. The HXP100 design also includes two programmable pivot points, enabling more flexibility to align a sample at a particular point or points of that sample.

To further ensure positioning performance, High Accuracy versions (HA) are available with guaranteed accuracy values. This enables the use of a Newport Hexapod in positioning applications, where position accuracy is required. In addition to accuracy along an axis, the Pitch and Yaw deviations during axial motion are also monitored and guaranteed. When the HA

Hexapod is used with RightPath™, this combination achieves positioning performance close to standard Newport stages.

Features

- Integrated 6-axis positioner
- Light, compact and low-profile
- No moving cables
- High stiffness (particular in z)
- No accumulation of motion errors
- Virtual center of rotation, set by software
- RightPath™ trajectory control

The HXP controllers have the capacity to drive up to two additional Single Axis stages while also providing advanced features including instrument grade I/O's, hardware based input triggers, event triggers, high-speed on-the-fly data acquisition, fast TCP/IP communication, and integrated TCL programming language for on-board processes. The HXP100 Series also takes advantage of low-runout, RightPath trajectory capability.

All these features improve accuracy and throughput, making the programmer's life much easier.




Applications

- Optics and satellite assembly and testing
- AED simulation
- Astronomy
- Biotechnology, surgery
- X-Ray diffraction
- Micromachining, micro-manipulation



The HXP100 design also includes two programmable pivot points, enabling more flexibility to align a sample at a particular point or points of that sample. The Newport Hexapod can not only relocate the pivot point, but through our advanced technology, the entire coordinate system can be relocated. In addition, two user-definable coordinate systems are


provided, called tool (moves with the Hexapod) and work (stationary coordinate systems). Incremental displacements are possible in either one of these systems in user-friendly Cartesian coordinates, and positions can be easily calculated from one system to the other by a function call, without the need for complex external coordinate transformations.



Absolute moves and positions are defined in the work coordinate system. Incremental moves can be done in the tool or in the work coordinate systems.



As a standard feature, the HXP controller allows the user to define a virtual pivot point in space for all rotations



RightPath™ Trajectory Control enables minimal runout in linear and arc trajectories.

Specifications

	HX100-MECA	HXP100HA-MECA	HXP100P-MECA	HXP100PHA-MECA	HXP100V6-MECA
Travel Range X, Y, Z ⁽¹⁾	±27.5, ±25, ±14 mm	±27.5, ±25, ±14 mm	±27.5, ±25, ±14 mm	±27.5, ±25, ±14 mm	±27.5, ±25, ±14 mm
Travel Range Θ X, Θ Y, Θ Z	±11.5, ±10.5, ±19 °	±11.5, ±10.5, ±19 °	±11.5, ±10.5, ±19 °	±11.5, ±10.5, ±19 °	±11.5, ±10.5, ±19 °
Minimum Incremental Motion X, Y, Z ⁽²⁾	0.5, 0.5, 0.25 μ m	0.50, 0.50, 0.25 μ m	0.10, 0.10, 0.05 μ m	0.10, 0.10, 0.05 μ m	0.5, 0.5, 0.25 μ m
Minimum Incremental Motion Θ X, Θ Y, Θ Z	0.25, 0.25, 0.5 mdeg	0.25, 0.25, 0.5 mdeg	0.05, 0.05, 0.10 mdeg	0.05, 0.05, 0.10 mdeg	0.25, 0.25, 0.5 mdeg
Uni-directional Repeatability X, Y, Z, Typical	±0.25, ±0.25, ±0.125 μ m	±0.14, ±0.13, ±0.05 μ m	±0.10, ±0.10, ±0.05 μ m	±0.10, ±0.10, ±0.05 μ m	±0.50, ±0.50, ±0.50 μ m
Uni-directional Repeatability X, Y, Z, Guaranteed	–	±0.25, ±0.25, ±0.125 μ m	–	±0.15, ±0.15, ±0.075 μ m	–
Uni-directional Repeatability Θ X, Θ Y, Θ Z, Typical	±0.125, ±0.125, ±0.25 mdeg	±0.125, ±0.125, ±0.25 mdeg	±0.05, ±0.05, ±0.10 mdeg	±0.05, ±0.05, ±0.10 mdeg	±0.25, ±0.25, ±0.50 mdeg
Accuracy XYZ, Guaranteed	–	±10, ±10, ±5 μ m	–	±5.0, ±5.0, ±2.5 μ m	–
Maximum Speed X, Y, Z	2.5, 2, 1 mm/s	2.5, 2, 1 mm/s	12, 10, 5 mm/s	12, 10, 5 mm/s	0.5, 0.5, 0.25 mm/s
Maximum Speed Θ X, Θ Y, Θ Z	1.8, 1.7, 3 °/s	1.8, 1.7, 3 °/s	8, 8, 16 °/s	8, 8, 16 °/s	0.2, 0.2, 0.4 °/s
Rigidity X, Y, Z ⁽³⁾	5, 5, 40 N/ μ m	5, 5, 40 N/ μ m	3, 3, 24 N/ μ m	3, 3, 24 N/ μ m	5, 5, 40 N/ μ m
Pitch X, Y, Z, Guaranteed	–	±75, ±75, ±75 μ rad	–	±37.5, ±37.5, ±37.5 μ rad	–
Yaw X, Y, Z, Guaranteed	–	±75, ±75, ±75 μ rad	–	±37, ±37, ±37 μ rad	–
Centered Load Capacity (4)	200 N	200 N	60 N	60 N	200 N
Cable Length	1.5 m	1.5m	3 m	3 m	1.5 m
Motor	DC Servo	DC Servo	DC Servo	DC Servo	Stepper motor
Weight	7.2 kg	7.2 kg	7.2 kg	7.2 kg	7.2 kg

¹⁾ Travel ranges are interdependent. The listed values are max. travels per axis when all other axes are in their centered position.


²⁾ Open loop values shown.

³⁾ Stiffness depends on Hexapod position. Values are given for all axes in their centered position.

⁴⁾ For Value shown for horizontal base plate. See graphs for maximum payload height and cantilever distance on next page

Max. Cantilever Distance of the Load



HXP100 Horizontal Base Plate



HXP100P Horizontal Base Plate



Dimensional Drawings




HXP100-MECA and HXP100HA-MECA


Note: Other top plate hole patterns or a center aperture are available upon request.




HXP100P-MECA Hexapod



Dimensional Drawings



HXP100P-MECA and HXP100P-MECA



Note: Other top plate hole patterns or a center aperture are available upon request.

Ordering Information

Model	Description
HXP100-MECA	Hexapod, 200 N load capacity
HXP100-ELEC-D ⁽¹⁾	Hexapod controller for HXP100-MECA
HXP100P-MECA	Hexapod Precision, 60 N load capacity
HXP100P-ELEC-D ⁽¹⁾	Hexapod controller for HXP100P-MECA

¹⁾ Contact Newport for the two additional SingleAxis drive capability

Model	Description
HXP100HA-MECA	Hexapod with guaranteed specifications, 200 N load capacity
HXP100HA-ELEC-D ⁽¹⁾	Hexapod controller for HXP100HA-MECA
HXP100PHA-MECA	Hexapod Precision with guaranteed specifications, 60 N load capacity
HXP100PHA-ELEC-D ⁽¹⁾	Hexapod controller for HXP100PHA-MECA

¹⁾ Contact Newport for the two additional SingleAxis drive capability

Note: Call Newport for quotes on the 10-6 hPa vacuum version.