



MFA Series

Miniature Linear Stages



 **Newport**[®]

USER'S MANUAL

Warranty

Newport Corporation warrants this product to be free from defects in material and workmanship for a period of 1 year from the date of shipment. If found to be defective during the warranty period, the product will either be repaired or replaced at Newport's discretion.

To exercise this warranty, write or call your local Newport representative, or contact Newport headquarters in Irvine, California. You will be given prompt assistance and return instructions. Send the instrument, transportation prepaid, to the indicated service facility. Repairs will be made and the instrument returned, transportation prepaid. Repaired products are warranted for the balance of the original warranty period, or at least 90 days.

Limitation of Warranty

This warranty does not apply to defects resulting from modification or misuse of any product or part.

CAUTION

Warranty does not apply to damages resulting from:

- **Incorrect usage:**
 - **Load on the stage greater than maximum specified load.**
 - **Carriage speed higher than specified speed.**
 - **Improper grounding.**
 - **Connectors must be properly secured.**
 - **When the load on the stage represents an electrical risk, it must be connected to ground.**
 - **Excessive or improper cantilever loads.**
- **Modification of the stage or any part thereof.**



CAUTION

Please return equipment in the original (or equivalent) packing.

You will be responsible for damage incurred from inadequate packaging if the original packaging is not used.

This warranty is in lieu of all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for a particular use. Newport Corporation shall not be liable for any indirect, special, or consequential damages.

© 2018 by Newport Corporation, Irvine, CA. All rights reserved.

Original instructions.

No part of this document may be reproduced or copied without the prior written approval of Newport Corporation. This document is provided for information only, and product specifications are subject to change without notice. Any change will be reflected in future publishings.

Table of Contents

Warranty	ii
EC Declaration of Conformity.....	v
Definitions and Symbols.....	vii
Warnings	viii
Caution	ix
<hr/>	
1.0 — Introduction.....	1
<hr/>	
2.0 — Description	2
2.1 Design Details	2
<hr/>	
3.0 — Characteristics.....	3
3.1 Definitions	3
3.2 Mechanical Specifications	4
3.3 Load Specification Definitions.....	4
3.4 Load Characteristics and Stiffness	5
3.5 Stage Weights	5
<hr/>	
4.0 — Drives and Motors	6
4.1 Stepper Drive Version	6
4.2 DC-Motor Drive Version	6
4.3 Sensor Position.....	7
4.4 Feedback Signal Position	8
4.5 Pinouts.....	8
<hr/>	
5.0 — Connection to Newport Controllers.....	10
5.1 Warnings on Controllers	10
5.2 Connection.....	11
5.3 Cables	11
<hr/>	
6.0 — Connection to Non-Newport Electronics	12
6.1 MFA-PPD & MFA-CC Connections	12
6.2 MFA-PP Stage	13
<hr/>	
7.0 — Dimensions.....	13
<hr/>	
8.0 — Accessories.....	14
8.1 MFA-BK Plate	14
8.2 (M-)MFA-TP Top Plates	14
8.3 MFA-BP Base Plate	14

9.0 — Maintenance	15
9.1 Maintenance	15
9.2 Repair	15
9.3 Calibration	15

Service Form	17
--------------------	----

EC Declaration of Conformity



2 Tech Drive
Andover, MA 01810
www.mksinst.com

EU27 Declaration of Conformity

Application of Council Directive(s):

- Electromagnetic Compatibility Directive (EMCD) – 2014/30/EU
- Machinery Directive – 2006/42/EC
- Restriction of Hazardous Substances Directive (RoHS2) – 2011/65/EU ⁽⁷⁾
- Restriction of Hazardous Substances Directive (RoHS3) – (EU) 2015/863⁽⁷⁾
- Waste Electrical and Electronic Equipment – Directive 2012/19/EU



Standard(s) to which conformity is declared:

- EN 61326-1:2013 – (EMC)
- EN ISO 12100:2010 Safety of Machinery – General Principles of Design – Risk Assessment and Risk Reduction

Emissions:

- EN 55011: 2016 +A1:2017 ⁽⁴⁾ Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement

Immunity:

- EN 61000-4-2:2009 EMC/Electrostatic Discharge Immunity Test
- EN 61000-4-3:2006+A2:2010 EMC/Radiated Radio Frequency Electromagnetic Field Immunity Test
- EN 61000-4-4:2012 EMC/Electrical Fast Transient/Burst Immunity Test
- EN 61000-4-5:2014+A1:2017 EMC/Surge Immunity Test
- EN 61000-4-6:2014 EMC/Conducted Disturbances induced by Radio Frequency Fields Immunity Test

Manufacturers Name: MKS Instruments, Inc., 2 Tech Drive, Andover, MA 01810 USA

Authorized Representatives Name & Location: _____

Equipment Type/Description: **Miniature Linear Stages**

Model Number(s) ⁽⁶⁾: **MFA-CC/-CCV6/-PP/-PPD.**

The object of the declaration described above is in conformity with the relevant Community harmonization legislation. MKS product conforms to the above Directive(s) and Standard(s) only when installed in accordance with manufacturer's specifications. This declaration has been issued under the sole responsibility of the manufacturer.

Date: 7/20/2021

Le Cointe Hervé – Quality Director

-
- 4) Class A, Group 2
 - 6) Compliance of the above model numbers requires the use of a braided shielded cable properly terminated at both ends – if so noted in the MKS Instruction Manual.
 - 7) RoHS Directive has to be checked for in scope products; cannot CE mark without compliance to RoHS. RoHS Directive can be unchecked only for systems which MKS sells which qualify for "Large Scale Industrial Tool" exclusion.

UK Declaration of Conformity



2 Tech Drive
Andover, MA 01810
www.mksinst.com

UK Declaration of Conformity

Application of Council Directive(s):

- Electromagnetic Compatibility Directive (EMCD) – 2014/30/EU
- Machinery Directive – 2006/42/EC
- Restriction of Hazardous Substances Directive (RoHS2) – 2011/65/EU ⁽⁷⁾
- Restriction of Hazardous Substances Directive (RoHS3) – (EU) 2015/863⁽⁷⁾
- Waste Electrical and Electronic Equipment – Directive 2012/19/EU



Standard(s) to which conformity is declared:

- BS EN 61326-1:2013 – (EMC)
- BS EN ISO 12100:2010 Safety of Machinery – General Principles of Design – Risk Assessment and Risk Reduction

Emissions:

- EN 55011: 2016 +A1:2017 ⁽⁴⁾ Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement

Immunity:

- BS EN 61000-4-2:2009 EMC/Electrostatic Discharge Immunity Test
- BS EN 61000-4-3:2006+A2:2010 EMC/Radiated Radio Frequency Electromagnetic Field Immunity Test
- BS EN 61000-4-4:2012 EMC/Electrical Fast Transient/Burst Immunity Test
- BS EN 61000-4-5:2014+A1:2017 EMC/Surge Immunity Test
- BS EN 61000-4-6:2014 EMC/Conducted Disturbances induced by Radio Frequency Fields Immunity Test

Manufacturers Name: MKS Instruments, Inc., 2 Tech Drive, Andover, MA 01810 USA

Authorized Representatives Name & Location: _____

Equipment Type/Description: **Miniature Linear Stages**

Model Number(s) ⁽⁶⁾: **MFA-CC/-CCV6/-PP/-PPD.**

The object of the declaration described above is in conformity with the relevant Community harmonization legislation. MKS product conforms to the above Directive(s) and Standard(s) only when installed in accordance with manufacturer's specifications. This declaration has been issued under the sole responsibility of the manufacturer.

Date: 7/20/2021

Le Cointe Hervé – Quality Director

4) Class A, Group 2

6) Compliance of the above model numbers requires the use of a braided shielded cable properly terminated at both ends – if so noted in the MKS Instruction Manual.

7) RoHS Directive has to be checked for in scope products; cannot CE mark without compliance to RoHS. RoHS Directive can be unchecked only for systems which MKS sells which qualify for "Large Scale Industrial Tool" exclusion.

Definitions and Symbols

The following terms and symbols are used in this documentation and also appear on the product where safety-related issues occur.

General Warning or Caution



The exclamation symbol may appear in warning and caution tables in this document. This symbol designates an area where personal injury or damage to the equipment is possible.

The following are definitions of the Warnings, Cautions and Notes that may be used in this manual to call attention to important information regarding personal safety, safety and preservation of the equipment, or important tips.



WARNING

Warning indicates a potentially dangerous situation which can result in bodily harm or death.



CAUTION

Caution indicates a potentially hazardous situation which can result in damage to product or equipment.

NOTE

Note indicates additional information that must be considered by the user or operator.

European Union CE Mark



The presence of the CE Mark on Newport Corporation equipment means that it has been designed, tested and certified as complying with all applicable European Union (CE) regulations and recommendations.

Warnings and Cautions



ATTENTION

This stage is a Class A device. In a residential environment, this device can cause electromagnetic interference. In this case, suitable measures must be taken by the user.

Warnings



WARNING

The motion of objects of all types carries potential risks for operators. Ensure the protection of operators by prohibiting access to the dangerous area and by informing the personnel of the potential risks involved.

WARNING

Do not use this stage when its motor is emitting smoke or is unusually hot to the touch or is emitting any unusual odor or noise or is in any other abnormal state.

Stop using the stage immediately, switch off the motor power and then disconnect the electronics power supply.

After checking that smoke is no longer being emitted contact your Newport service facility and request repairs. Never attempt to repair the stage yourself as this can be dangerous.

WARNING

Make sure that this stage is not exposed to moisture and that liquid does not get into the stage.

Nevertheless, if any liquid has entered the stage, switch off the motor power and then disconnect the electronics from power supply.

Contact your Newport service facility and request repairs.



WARNING

Do not insert or drop objects into this stage, this may cause an electric shock, or lock the drive.

Do not use this stage if any foreign objects have entered the stage. Switch off the motor power and then disconnect the electronics power supply.

Contact your Newport service facility for repairs.

WARNING

Do not place this stage in unstable locations such as on a wobbly table or sloping surface, where it may fall or tip over and cause injury.

If this stage has been dropped or the case has been damaged, switch off the motor power and then disconnect the electronics power supply.

Contact your Newport service facility and request repairs.

WARNING

Do not attempt to modify this stage; this may cause an electric shock or downgrade its performance.

WARNING

Do not exceed the usable depth indicated on the mounting holes (see section “Dimensions”). Longer screws can damage the mechanics or cause a short-circuit.

Caution

CAUTION

Do not place this stage in a hostile environment such as X-Rays, hard UV,... or in any vacuum environment.

Only the MFA-CCV6 stage is compatible and can be used in a vacuum environment up to 10^{-6} hPa.

CAUTION

Do not place this stage in a location affected by dust, oil fumes, steam or high humidity. This may cause an electric shock.

CAUTION

Do not leave this stage in places subject to extremely high temperatures or low temperatures. This may cause an electric shock.

- Operating temperature: +10 to +35 °C
- Storage temperature: -10 to +40 °C (in its original packaging)

CAUTION

Do not move this stage if its motor power is on.

Make sure that the cable to the electronics is disconnected before moving the stage. Failure to do so may damage the cable and cause an electrical shock.

CAUTION

Be careful that the stage is not bumped when it is being carried. This may cause it to malfunction.

CAUTION

When handling this stage, always unplug the equipment from the power source for safety.

CAUTION

When the carriage is in its end-of-run position, it is strongly recommended not to go beyond this point as this may damage the stage mechanism.

CAUTION

Contact your Newport service facility to request cleaning and specification control every year.

CAUTION

Do not manually push or pull the carriage. You may easily exceed acceptable axial load and damage the stage.



Miniature Linear Stages

MFA Series

1.0 Introduction

This manual provides operating instructions for the stage that you have purchased in the MFA Series:

- MFA-CC
- MFA-CCV6 ⁽¹⁾
- MFA-PP
- MFA-PPD

¹⁾ REMARK

Vacuum compatible stage to 10^{-6} hPa. In this case, the max. speed and load capacity have to be divided by two.



MFA-CC miniature linear stage.

RECOMMENDATION

We recommend you read carefully the chapter “Connection to electronics” before using the MFA stage.



MFA stages in an XYZ configuration.

2.0 Description

Designed for space-limited applications and compact multi-axis assemblies, MFA Series linear stages supply very high resolution, single-axis translation in a low-cost, motorized package. Typical applications for this stage are fiber optic alignment, laser diode research, bio-medical applications and inspection systems.

MFA linear stages are available in two versions: The MFA-CC with a DC-motor, features an integrated gear-box and a motor mounted high resolution 2,048 cts/rev rotary encoder. The high-resolution position feedback and low-friction mechanical design ensures ultra-smooth motion with 100 nm sensitivity. The DC motor supplies an optimized output torque that permits the use of a lower ratio step-down gear allowing for faster motion with higher reliability and lower backlash. Hence, the MFA-CC is the recommended choice for applications that require small incremental motion with high dynamic speed range and good repeatability. The MFA-PP and MFA-PPD stepper motor versions are more economical solutions for less demanding applications.

Travel limit switches prevent bearing damage from accidental over-travel.

MFA stages feature an all-steel construction that provides a higher stiffness-to-weight ratio and lower thermal expansion compared to aluminum designs. The result is superior performance in a smaller footprint. The smooth motion of the MFA linear stages is further accentuated by Newport's proprietary double-row linear ball bearing design with bearing ways that are directly machined into the structural frame of the stage. Compared to other solutions that use commercial bearings, MFA linear stages have a higher load capacity and stiffness with low pitch and yaw errors.

Another benefit of Newport's integrated bearing ways is the availability of 4 widely spaced mounting holes for base mounting and XY-assemblies. This provides better stress distribution and allows for more rigid multi-axis combinations than other designs that provide only line contact with 2 mounting holes.

2.1 Design Details

Base Material	Stainless Steel
Bearings	Double row linear ball bearings
Drive Mechanism	Backlash compensated leadscrew
Drive Screw Pitch	0.5 mm
Reduction Gear	MFA-CC: 1:14 MFA-PP: 1:43
Feedback	MFA-CC: Motor mounted rotary encoder; 2,048 cts/rev MFA-PP: None
Limit Switches	Optical switches
Origin	Uses motor side limit for homing, typically $\pm 2.5 \mu\text{m}$ repeatability
Cable Length	MFA-PP, MFA-PPD & MFA-CC: 3 m MFA-CCV6: 1.5 m

3.0 Characteristics

3.1 Definitions

Specifications of our products are established in reference to ISO 230 standard part II “Determination of accuracy and repeatability of positioning numerically controlled axes”.

This standard gives the definition of position uncertainty which depends on the 3 following parameters:

Absolute Accuracy

Difference between ideal position and real position.

Accuracy

Difference between ideal position and real position after the compensation of linear errors.

Linear errors include: cosine errors, inaccuracy of screw or linear scale pitch, angular deviation at the measuring point (Abbe error) and thermal expansion effects. All Newport motion electronics can compensate for linear errors.

The relation between absolute accuracy and on-axis accuracy is as follows:

$$\text{Absolute Accuracy} = \text{Accuracy} + \text{Correction Factor} \times \text{Travel}$$

Repeatability

Ability of a system to achieve a commanded position over many attempts.

Reversal Value (Hysteresis)

Difference between actual position values obtained for a given target position when approached from opposite directions.

Minimum Incremental Motion (MIM or Sensitivity)

The smallest increment of motion a device is capable of delivering consistently and reliably.

Resolution

The smallest increment that a motion device can theoretically move and/or detect. Resolution is not achievable, whereas MIM, is the real output of a motion system.

Yaw, Pitch

Rotation of carriage around the Z axis (Yaw) or Y axis (Pitch), when it moves.

The testing of accuracy, repeatability, and reversal error are made systematically with test equipment in controlled environment (20^{±1} °C).

A linear cycle with 21 data points on the travel and 4 cycles in each direction gives a total of 168 points.

Guaranteed and Typical Specifications

Guaranteed maximum performance values are verified per Newport's A167 metrology test procedure. For more information, please consult the metrology tutorial section in the Newport catalog or at www.newport.com

3.2 Mechanical Specifications



	MFA-PP and MFA-PPD	MFA-CC
Travel Range [in. (mm)]		1 (25)
Minimum Incremental Motion (μm)	0.1	0.1
Uni-directional Repeatability, Typical (Guaranteed) (μm)	0.12 (±0.25)	±0.08 (±0.15)
Bi-directional Repeatability, Typical (Guaranteed) ⁽¹⁾⁽²⁾ (μm)	±0.2 (±0.75)	±0.15 (±0.75)
Accuracy, Typical (Guaranteed) ⁽¹⁾ (μm)	±0.9 (±3.0)	±0.7 (±3.0)
Maximum Speed (mm/s)	0.3 (MFA-PP) 1.0 (MFA-PPD)	2.5
Pitch, Typical (Guaranteed) ⁽¹⁾ (μrad) ⁽³⁾		±25 (±100)
Yaw, Typical (Guaranteed) ⁽¹⁾ (μrad) ⁽³⁾		±30 (±50)
MTBF	10,000 h at a 1 kg load with a 20% duty cycle	

¹⁾ For the definition of Typical and Guaranteed specifications see "Motion Basics Terminology & Standards" Tutorial at www.newport.com

²⁾ After backlash compensation.

³⁾ To obtain arcsec units, divide μrad value by 4.8.

NOTE

For the MFA-CCV6 vacuum compatible stage to 10⁶ hPa, the max. speed and load capacity have to be divided by two.



CAUTION

To reach specifications stated, stages must be fixed on a plane surface with a flatness of 5 μm.

3.3 Load Specification Definitions

Normal Load Capacity (Cz)

Maximum load a stage can move while maintaining specifications.

This value is given with speed and acceleration specified for each stage, and with a load perpendicular to bearings.

Axial Load Capacity (±Cx)

Maximum load along the direction of the drive train.


Off-Centered Load (Q)

Maximum cantilever-load a stage can move: $Q \leq Cz / (1 + D/a)$

D: Cantilever distance.

a: Construction parameter

3.4 Load Characteristics and Stiffness



C_z , Normal centered load capacity (N)	50
$-C_x, +C_x$, Axial load capacity (N)	10
$K_{\alpha x}$, Compliance in roll ($\mu\text{rad}/\text{Nm}$)	60
$K_{\alpha y}$, Compliance in pitch ($\mu\text{rad}/\text{Nm}$)	10
Q , Off-center load (N)	$Q \leq C_z \div (1 + D/20)$
Where D = Cantilever distance (mm)	

NOTE

For the MFA-CCV6 vacuum compatible stage to 10^{-6} hPa, the load capacity has to be divided by two.

3.5 Stage Weights

The weight indicated is the average value for stages with their cable installed.

	Weight [lb (kg)]
MFA	1.3 (0.6)

The weight variation between drive units is not very significant.



MFA-CCV6 stage.

4.0 Drives and Motors

4.1 Stepper Drive Version

Stepper-motor-driven stages are offered in one half-step drive version: MFA-PP & MFA-PPD.

Micro-Step Drive


This is the drive for stepper or pulse-driven motors, transmitted by the electronic unit, which entails a theoretical movement of the motor by one fraction of a full-step. For these translation stages, the micro-step is equivalent to 1/64 of the full-step.

Stepper Motor Performance Specifications and Characteristics


	Resolution ⁽¹⁾ (μm)	Speed (mm/s)	Angle by Step ($^\circ$)	RMS Current per Phase (A)	Resistance (Ω)	Inductance (mH)
MFA-PP	0.00757	0.3	15	0.25	12.5	5.5
MFA-PPD	0.00757	1.0	15	0.25	12.5	5.5

¹⁾ Micro-step driving.

Command Signals for the Stepper Motor



Step	1	2	3	4	5	8	7	8
+ Phase 1	+	+	-	-	+	+	-	-
- Phase 1	-	-	+	+	-	-	+	+
+ Phase 2	-	+	+	-	-	+	+	-
- Phase 2	+	-	-	+	+	-	-	+
Direction	-	←	MOTION	→	+			



4.2 DC-Motor Drive Version

One DC-motor-driven configuration is available: MFA-CC.


DC-Motor Performance Specifications and Characteristics

	Resolution (μm)	Speed (mm/s)	Nominal Voltage (V)	Max RMS Current (A)	Max. Peak Current (A)	Resistance (Ω)	Inductance (mH)	Tachometer Const. (V/krpm)
MFA-CC	0.0174	2.5	24	0.15	0.21	54.6	1.19	-

NOTE

For the MFA-CCV6 vacuum compatible stage to 10^{-6} hPa, the max. speed has to be divided by two.


Command Signals for the DC-Motor



In the above drawings, + Motor signal is referred to - Motor signal.

- ① When the stage moves in + Direction, the + Motor voltage is higher than - Motor voltage.
- ② When the stage moves in - Direction, the + Motor voltage is lower than - Motor voltage.

4.3 Sensor Position




End-of-Run are 5 V open collector type.




CAUTION

“End-of-Run” are active signals and should not be connected to any other source.

4.4 Feedback Signal Position




The incremental sensor consists of an optical scale and an encoder head. When the carriage moves, the encoder head generates square signals in quadrature and sends to pins #19, #20, #23 and #24 of the SUB-D25 connector.



“Encoder” are “differential pair” (type RS-422) type output signals. Using these signals permits a high immunity to noise. Emission circuits generally used by Newport are 26LS31 or MC3487. Reception circuits to use are 26LS32 or MC3486.


4.5 Pinouts

The pinout diagram for the MFA-PP stages mini-DIN9 connector is shown below.



MFA-PP	
1	- End-of-Run
2	Phase 1a
3	+5 V
4	Phase 2b
5	N.C.
6	Phase 2a
7	Phase 1b
8	Ground
9	+ End-of-Run

The pinout diagrams for the MFA-PPD and MFA-CC stages SUB-D25M connector are shown below.



MFA-PPD		MFA-CC	
1	Phase 1a	1	N.C.
2	N.C.	2	N.C.
3	Phase 1b	3	N.C.
4	N.C.	4	N.C.
5	Phase 2a	5	+ Motor
6	N.C.	6	+ Motor
7	Phase 2b	7	- Motor
8	N.C.	8	- Motor
9	N.C.	9	N.C.
10	N.C.	10	N.C.
11	N.C.	11	N.C.
12	N.C.	12	N.C.
13	Reserved ⁽¹⁾	13	Reserved ⁽¹⁾
14	Ground	14	Ground
15	N.C.	15	N.C.
16	0 V	16	0 V
17	+ End-of-Run	17	+ End-of-Run
18	- End-of-Run	18	- End-of-Run
19	N.C.	19	Encoder Phase A
20	N.C.	20	Encoder Phase B
21	+5 V	21	+5 V
22	0 V	22	0 V
23	N.C.	23	Encoder Phase /A
24	N.C.	24	Encoder Phase /B
25	N.C.	25	N.C.

¹⁾ Pin #13 of SUB-D25M is reserved for internal logic for Newport controller.

5.0 Connection to Newport Controllers

5.1 Warnings on Controllers

Controllers are intended for use by qualified personnel who recognize shock hazards and are familiar with safety precautions required to avoid possible injury. Read the controller user's manual carefully before operating the instrument and pay attention to all written warnings and cautions.

WARNING

Disconnect the power plug under the following circumstances:

- If the power cord or any attached cables are frayed or damaged in any way.
- If the power plug is damaged in any way.
- If the unit is exposed to rain, excessive moisture, or liquids are spilled on the unit.
- If the unit has been dropped or the case is damaged.
- If you suspect service or repair is required.
- Whenever you clean the electronics unit.

CAUTION

To protect the unit from damage, be sure to:

- Keep all air vents free of dirt and dust.
- Keep all liquids away from the unit.
- Do not expose the unit to excessive moisture (85% humidity).
- Read this manual before using the unit for the first time.



WARNING

All attachment plug receptacles in the vicinity of this unit are to be of the grounding type and properly polarized.

Contact your electrician to check your receptacles.

WARNING

This product is equipped with a 3-wire grounding type plug.

Any interruption of the grounding connection can create an electric shock hazard.

If you are unable to insert the plug into your wall plug receptacle, contact your electrician to perform the necessary alterations to ensure that the green (green-yellow) wire is attached to earth ground.

WARNING

This product operates with voltages that can be lethal.

Pushing objects of any kind into cabinet slots or holes, or spilling any liquid on the product, may touch hazardous voltage points or short out parts.

5.2 Connection

There is a label on every stage indicating its part and serial numbers.



WARNING

Always turn the controller's power OFF before connecting to a stage.

NOTE




These stages are ESP compatible. Enhanced System Performance is Newport's exclusive technology that enables Newport ESP motion controllers to recognize the connected Newport ESP stage and upload the stage parameters. This ensures that the user can operate the motion system quickly and safely.

5.3 Cables

MFA-PP, MFA-PPD and MFA-CC are delivered equipped with a 3-meter cable, the MFA-CCV6 stage is delivered with a 1.5-meter cable.

These cables are equipped with a SUB-D25M connector for direct connection to Newport Controllers.



WARNING

This cable is shielded correctly. For a correct operation, make sure to lock connectors (ground continuity provided by the cable).



WARNING

Keep the motor cables at a safe distance from other electrical cables in your environment to avoid potential cross talk.

6.0 Connection to Non-Newport Electronics

6.1 MFA-PPD & MFA-CC Connections

WARNING

Newport is not responsible for malfunction or damage of MFA stages when used with non-Newport controllers.

WARNING


Newport guarantees “ C ” compliance of MFA stages only if used with Newport cables and controllers.

It is the customer’s responsibility to modify the cable and take care of sensor signal connections, when using the stage with non-Newport controllers.




End-of-Runs and Mechanical Zero are open collector type with a 5.6 V protective Zener diode.

For connection to non-Newport electronics, pin #17 and #18 End-of-Runs must be connected to +5 V with a pull-up resistor. For some non-Newport controllers, the +5 V output signal may be pulled up internally.



To enable End-of-Run, when a MFA stage is used with non-Newport controllers using the Sub-D25M connector provided, pin #13 must be connected to +5 V with a pull-up resistor (ex: 4.99 k Ω).



6.2 MFA-PP Stage




WARNING

The MFA-PP stage is equipped with a mini-DIN9 non-standard connector and must only be used with our Newport NSC200 controller/driver.



NSC200 Single-Axis Low-Power Motion Controller/Driver.

7.0 Dimensions




8.0 Accessories

NOTE


These plates are not included with MFA stages and must be ordered separately.

8.1 MFA-BK Plate


The MFA-BK top plate is used for XZ and XYZ assemblies with MFA stages.



8.2 (M-)MFA-TP Top Plates



8.3 MFA-BP Base Plate



9.0 Maintenance

RECOMMENDATION

Please contact Technical Sales Support team for recommendations on application specific maintenance.

9.1 Maintenance

The MFA stage requires no particular maintenance. Nevertheless, this is a precision mechanical device that must be kept and operated with caution.

PRECAUTIONS

The MFA stage must be used or stocked in a clean environment, without dust, humidity, solvents or other substances.

RECOMMENDATION

It is recommended to return the stage to Newport for re-lubrication after 2000 hours of use.

If the MFA stage is mounted on a workstation and cannot be easily removed, please contact Newport's After Sales Service for further instructions.

9.2 Repair

CAUTION



Never attempt to disassemble a component of the stage that has not been covered in this manual.

To disassemble a non specified component can cause a malfunction of the stage.

If you observe a malfunction in your stage, please contact us immediately to arrange for a repair.



CAUTION

Any attempt to disassemble or repair a stage without prior authorization will void your warranty.

9.3 Calibration



CAUTION

It is recommended to return your MFA stage to Newport once a year for recalibration to its original specifications.



Visit MKS | Newport Online at:
www.newport.com

North America & Asia

Newport Corporation
1791 Deere Ave.
Irvine, CA 92606, USA

Sales

Tel.: +1 (949)-863-3144
e-mail: sales@newport.com

Technical Support

Tel.: +1 (949)-863-3144
e-mail: tech@newport.com

Service, RMAs & Returns

Tel.: +1 (949)-863-3144
e-mail: service@newport.com

Europe

MICRO-CONTROLE Spectra-Physics S.A.S
7 rue des Plantes
45340 Beaune-la-Rolande
France

Sales Europe (EMEA)

Tel.: +49 (0) 6151-708-0
e-mail: germany@newport.com

Sales France

Tel.: +33 (0)1 60 91 68 68
e-mail: france@newport.com

Sales UK

Tel.: +44 (0)1235 432 710
e-mail: uk@newport.com

Technical Support

e-mail: tech_europe@newport.com

Service & Returns

Tel.: +33 (0)2 38 40 51 55
DST-BEA-RMA-service@newport.com