



CONEX-LDS

Electronic Autocollimator



USB Driver Installation Manual

©2017 by Newport Corporation, Irvine, CA. All rights reserved.

Original instructions.

No part of this document may be reproduced or copied without the prior written approval of Newport Corporation. This document is provided for information only, and product specifications are subject to change without notice. Any change will be reflected in future publishings.

Table of Contents

1.0 Overview and Setup.....	1
1.1 Components	1
1.2 Electrical Installation	3

2.0 First Connection	4
2.1 CONEX-LDS USB Driver Installation on Windows 7.....	4
2.2 USB Device Installation for Windows XP	5

3.0 CONEX-LDS Installation Verification.....	7
---	----------

Service Form	9
---------------------------	----------



Electronic Autocollimator CONEX-LDS

1.0 Overview and Setup

1.1 Components

USB Cable



CONEX-USB-RS422 – Serial to USB Interface Converter



CONEX-LDS-CABLEXX



CONEX-LDS



CONEX-LDS-PS



1.2 Electrical Installation

1. Connect the LDS cable and power supply cable to the LDS head
2. Connect the LDS cable and USB cable to the serial converter



3. Connect the USB cable to an available USB port of the PC



2.0 First Connection




CAUTION

BEFORE ANY INSTALLATION, CONTACT YOUR I.T. ADMINISTRATOR TO VERIFY THAT YOU HAVE THE APPROPRIATE RIGHTS.


2.1 CONEX-LDS USB Driver Installation on Windows 7

1. Connect the CONEX-LDS controller to a USB port with the provided USB cable
2. Detection of this new connected device (the first time) is signaled by a message in the bottom -right of the screen




3. For Windows 7 64 bit, the USB driver installation is completed when the message disappears
4. For Windows 7 32 bits, right click on the new detected device and select “Update driver”.

Click “Search automatically for updated driver software”



5. The following window will pop up after driver installation completes



2.2 USB Device Installation for Windows XP

1. Download CONEX-USB-RS422_driver.zip from newport.com
2. Run USBCom-CDM_20828.exe



3. Select “Extract” and to following window will pop up



The following page indicates that the drivers have been successfully installed



4. Once the drivers are installed connect the CONEX-LDS to a USB port with the provided USB cable
5. Detection of this new connected device (the first time) is signaled by a message in the bottom -right of the screen




Now, your CONEX-LDS is installed and ready to use


3.0 CONEX-LDS Installation Verification

Open the “Device Manager” and select the “Device” tab to check your configuration. In the display menu, select “View devices by type”

Expand “USB Bus controllers” line and check the new line USB Serial Converter:



Expand “Ports (COM and LPT)” line and check the new line USB Serial Port (COMx):



NOTES

Note the COM port number (x) assigned to the device. This will be required for troubleshooting.



Visit Newport Online at:
www.newport.com

North America & Asia

Newport Corporation
1791 Deere Ave.
Irvine, CA 92606, USA

Sales

Tel.: (800) 222-6440
e-mail: sales@newport.com

Technical Support

Tel.: (800) 222-6440
e-mail: tech@newport.com

Service, RMAs & Returns

Tel.: (800) 222-6440
e-mail: service@newport.com

Europe

MICRO-CONTROLE Spectra-Physics S.A.S
9, rue du Bois Sauvage
91055 Évry CEDEX
France

Sales

Tel.: +33 (0)1.60.91.68.68
e-mail: france@newport.com

Technical Support

e-mail: tech_europe@newport.com

Service & Returns

Tel.: +33 (0)2.38.40.51.55

