

# **Model 8742 Picomotor™ Controller/Driver Firmware**

Version 3.04

Revision Date: September 09, 2016

## **New Focus World Wide Web Server**

You can access a variety of Newport Corporation information sources via the Newport Web Server, <http://www.newport.com>.

## **Technical Support**

Please review the User's Manual first if you are experiencing difficulties with the product. The following Technical Support information is listed below if you still need help.

Tel: 1-800-222-6440  
Fax: 1-949-253-1479  
Email: [rma.service@newport.com](mailto:rma.service@newport.com)  
Internet: Navigate for "Support" at <http://www.newport.com>

# What's New

The following is a list of added features and corrections to Model 8742 Picomotor™ Controller/Driver firmware:

## VERSION 3.04 September 09, 2016

*Fixes:* Allows a \n or \r character to terminate a command on the telnet interface, rather than just a \r character.

*Feature Enhancements:* Initial ship release version for TIVA-version product

## VERSION 3.03 September 02, 2016

*Fixes:* None

*Feature Enhancements:* Initial ship release version for TIVA-version product (retracted)

©2016 New Focus. All rights reserved.

This document is for informational purposes only. New Focus makes no warranties, express or implied, in this document.  
Microsoft and Windows are registered trademarks of Microsoft Corporation.  
IBM is a registered trademark of International Business Machines Corporation.