



PART NUMBER	A in [mm]
VIS-LSS-24	24 [600]
VIS-LSS-30	30 [750]
VIS-LSS-36	36 [900]



NOTES:

1. VIS-LSS-XX CAN BE ADJUSTED FOR HEIGHTS CORRESPONDING TO BREADBOARDS OF 2" THK AND 4" THK.
2. WORKSTATION IS SHOWN FOR DESMONSTRATION PURPOSE AND DOESN'T COME WITH VIS-LSS SHELF PRODUCT.
3. DIMS IN [] ARE IN mm

