



MODEL NUMBER	DESCRIPTION	LEVELING ARM	NO. OF BARB FITTING	FITTING CONNECTION
	ISOLATOR, STANDARD	YES	1	TO AIR SUPPLY MANIFOLD
CM-225-1S	ISOLATOR, SLAVE	NO	1	TO ISOLATOR, VALVED
CM-225-1V	ISOLATOR, MASTER	YES	2	ONE TO ISOLATOR, SLAVE ONE TO AIR SUPPLY MANIFOLD

MODEL NUMBER	DESCRIPTION
CM-225	SET OF FOUR ISOLATORS; TWO STANDARD UNITS, ONE MASTER UNIT, ONE SLAVE UNIT
CM-225-3	SET OF THREE STANDARD ISOLATORS
UCM-225	SET OF FOUR CLEANROOM ISOLATORS, TWO STANDARD UNITS, ONE MASTER UNIT, ONE SLAVE UNIT

NOTES:

1. A STANDARD UNIT HAS ITS OWN VALVE TO CONTROL FLOATING.
2. A SLAVE UNIT IS PNEUMATICALLY HOOKED UP TO A MASTER UNIT, WHICH CONTROLS FLOATING FOR BOTH UNITS.
3. DIMS IN [] ARE IN mm

